



3DMed Sample Transport Medium


(Extraction Free)



3DMed

Innovative Sample Preparation Solution

PCR detection with 3DMed Sample Transport Medium(Extraction Free)



Excellent Product Performance

One step operation

Inactivation together with lysis without adding nucleic acid releasing agent

Time saving reduces labor

The lysis products can be directly used in qPCR



Good stability

Stored steadily at room temperature for 1 day after sampling

High sensitivity

No washing and elution
No loss of nucleic acid

Product Specifications

Applicable Specimens	Swab
Tube volume	10ml/1.5ml
Liquid loading volume	1ml/0.5ml
Packaging	50 PCS/ box
Accessory	NP swabs included



Website

3DMed

Building 2A, 158 XinJunhuan Road,
Shanghai, China

(+86) 400-021-1661

ivdsupport@3dmedcare.com